

## CoproTi-5

**CE 0483**

### Specifications

#### **product:**

product type:

product shape:

CE-mark:

applied standards:

veneer porcelain:

**printing date: 13.02.2012**

#### **CoproTi-5**

Titan Grade 5 – milling blank

metal disc 98,3mm Ø in different diameters and thicknesses

CE 0483

DIN EN ISO 9001:2008, DIN ISO 5832-3 and ASTM F67  
manufacturing and testing according to DIN EN ISO 13485 and medical  
products guideline 93/42/EEC annex II excluding section 4  
all standard veneering porcelains for titanium

### **composition:**

**Ti** ~90

**Fe** 0,13

**Al** 6,0

**V** 4,1

**O** 0,1

**H** 0,001

**N** 0,01

**C** 0,03

### **mechanical properties:**

density: 4,45 g/cm<sup>3</sup>

vickers hardness: 353

module of elasticity: 114 kN/mm<sup>2</sup>

CTE: 836 MPa

yield strength: 897 MPa

tensile strength: 12 %

fracture strain: 10,3 10<sup>-6</sup>/m. Celsius

### description

CoproTi-5 is a high-quality titanium alloy (grade 5-ELI) for the CAD / CAM technology. This industrially manufactured material ensures consistent quality, has high tensile and hardness values

CoproTi-5 naturally is biocompatible.

### indication:

**CoproTi-2** (pure titanium Grade 2) single crowns, small bridges in anterior and posterior region

**CoproTi-4** (pure titanium Grade 4) single crowns, medium bridges, implants in anterior and posterior region, implants and bar constructions

**CoproTi-5** (titanium alloy Grade 5-ELI) single crowns up to big bridges and bar constructions in anterior and posterior region, implants

### Instruction for use:

- Cut out, smoothen frameworks and single elements with suitable milling burs for titanium.

### Cleaning:

- Fettle and smoothen the surfaces of milled frameworks with special, titanium suitable cross-cut burs or separating discs in only one direction to avoid a blistering in the porcelain
- Sandblast the frameworks with 110µ (2-3 bar pressure) aluminum oxide and steam clean or dip them in methylalcohol. Never use hydrofluoric acid!

### Bonding of ceramic:

- Remove oxides after firing by blasting with glass beads. Finish with rubber stones and polishing paste
- Please follow the instructions for use of your chosen veneering porcelain manufacturer

**Hazard note! During dry milling of titanium, chips and swarfs can ignite themselves and cause fire.**