

## CopranColor/ CopranUltra-T Color

CopranColor is a colouring liquid in 16 classic shades, A1-D4

CopranUltra-T Color is especially adapted to the characteristics of highly translucent zirconium

### your advantages at a glance

- water based, does not contain acid
- therefore does not corrode metal instruments or metal workbench surfaces
- does not corrode the heating elements in the sintering furnace
- does not contain oily substances. Therefore it can fully benefit from the capillary attraction of zirconium and infiltrate quickly and deeply into the structure
- cutting sinter drops, connecting bars or milling telescopes will not lead to lighter areas or blotches
- can easily be dissolved with distilled water to lighten the tone if the colours are too intensive for your taste
- can be used with all types of yttrium stabilised zirconium
- reproduces the 4 classic colour groups precisely

