



Copran Zr & Copran Zri

perfection for a lifetime – zirconium with a lifetime warranty*

Our zirconium blanks are exclusively made from ready to use, extra fine, spray granulated raw material from Tosoh / Japan – the world leader for technical zirconium. Due to its quality and our production in Germany, we are proud to offer a lifelong warranty for our zirconium products. Our high strength Copran Zr blanks are axial pressed, Copran Zri blanks are **isostatic manufactured**.

Copran Zr Ultra-T & Zri Ultra-T

is a high translucent material without the problem of low temperature degradation. Due to its contents of aluminium oxide it is long term stable and we also offer **our lifetime warranty** on that material as well.

composition

zirconium oxide (ZrO ₂):	95
yttrium oxide (Y ₂ O ₃):	5
aluminium oxide (Al ₂ O ₃):	<0,4 (Copran Zr & Copran Zri)
4-point bending strength:	MPa (N/mm ⁻²) 1327 to DIN EN ISO 6872
aluminium oxide (Al ₂ O ₃):	<0,2 (Copran Zr Ultra-T & Zri Ultra-T)

*terms of warranty at www.white-peaks-dental.com

